

ELISA Solutions

Premium Seafood Testing



Prawn/Shrimp



Fish



Basa

Drug Residue	Cat. No	Matrices	LOD
AHD FAST*	NF3463	 	0.08 ppb
AMOZ FAST*	NF3462	 	0.08 ppb
AOZ FAST*	NF3465	 	0.06 ppb
SEM FAST*	NF3461	 	0.25 ppb
Chloramphenicol	CN1469	 	0.02 ppb
Chloramphenicol FAST*	CN10171	 	0.02 ppb
Doxycycline Sensitive*	DXS10119 A/B	 	14 ppb
Enrofloxacin/Ciprofloxacin*	EC3473	 	Method 1: ≤ 0.75 ppb Method 2: ≤ 0.5 ppb
Malachite Green FAST*	MGI0108 A/B	 	Method 1: 0.5 ppb (prawn, shrimp, basa) Method 2: 0.25 ppb (prawn, shrimp)
DNSH*	DNS10592	 	0.25 ppb

*These kits provide analysis under 60 mins

ELISA Plate Reader

We recommend the following ELISA plate reader, offering optimal performance with every test.

800TS Plate Reader

The 800TS is a robust microplate reader ideally suited for routine drug detection in the laboratory. Complete with touchscreen and USB capabilities, the high-quality microplate reader requires limited space for outstanding performance, with Gen5™ software available for advanced data handling and analysis.

- Compact system
- Speed reading
- Optimal performance
- High accuracy
- Multiple plate formats
- Gen5™ software

