

INFINIPILEX

ENSURE MILK QUALITY | EVERY DROP | EVERY TIME

The Randox InfiniPlex range has been developed specifically for the dairy industry as a surveillance tool to improve and monitor milk quality and safety worldwide.

Ideal for Use In:



Dairy Processors



Government Food Testing Laboratories



Food Research Institutes

IPM I & II	IPM III	IPM IV and V
<p>Precise screening for thorough milk surveillance and adherence to EU regulatory standards.</p> <p>Accurate and fast detection of multiple antibiotic residues, contaminants and antiparasitics in real-time.</p>	<p>Detection of multiple prevalent pesticides.</p> <p>Key for ensuring milk quality, safeguarding health, meeting legal standards, ensuring consumer trust, and promoting sustainable agriculture.</p>	<p>Differentiated sensitivities for a more tailored screen, ideal for first-stage quality checks.</p> <p>Fast screening and results allow for timely data-driven decision-making.</p> <p>Detection of multiple antibiotic residues, contaminants and antiparasitics.</p>

Why Choose Randox?



The World's First

InfiniPlex II is the ONLY milk screening test available globally for the detection of multiple antiparasitic contaminants*.



Meets EU Regulation

InfiniPlex has a >90% compliance rate with all currently regulated antibiotics within the EU from a single test.



Tailored Sensitivity

Different kit options are available for both increased and reduced sensitivity levels to suit your testing requirements.



Vast Test Menu

Simultaneously detect up to 130 antibiotics and contaminants from a single milk sample with InfiniPlex.



Identify Exact Veterinary Drugs

Results allow you to pinpoint the exact veterinary drug administered at an animal level.



High Throughput

Run up to 48 milk samples for multiple contaminants in under 2 hours.

* Covers Both EU Authorised and Unauthorised Antiparasitic Veterinary Medicines

Available On:



Evidence Investigator

- Semi-Automated Benchtop Analyser
- Qualitative Results
- Results for 48 Samples in Under 2 Hours



Evidence MultiSTAT

- Fully Automated Benchtop Analyser
- Qualitative Results
- Results for 1 Sample in Under 30 Minutes



randoxfood.com



info@randoxfood.com